



EAGLE
BIOSCIENCES

Nitric Oxide Synthase Ultrasensitive Colorimetric Assay

Catalog Number: NSU03-K01

96 Wells

For Research Use Only

v. 1.1 (20JUN23)

EAGLE BIOSCIENCES, INC.

20A NW Blvd., Suite 112, Nashua, NH 03063
Phone: 866-419-2019 Fax: 617-419-1110

WWW.EAGLEBIO.COM




INTRODUCTION

The traditional method for measuring nitric oxide synthase (NOS) activity is performed by radiochemical assay that measures the conversion of L-[³H]arginine to L-[³H]citrulline. This method is expensive and requires regulation of radioactive materials. The Ultrasensitive Colorimetric NOS Assay Kit is a low-cost novel assay that allows for the detection of NOS activity without the need for radioactivity. Our Ultrasensitive NOS Assay Kit employs a NADPH recycling system to allow NOS to operate linearly for hours as nitric oxide-derived nitrate and nitrite accumulate. NOS can be assayed spectrophotometrically by measuring the accumulation of its stable degradation products, nitrate and nitrite. The ratio of these two products in biological fluids, tissue culture media, etc. may vary substantially. Hence, for accurate assessment of the total nitric oxide generated, one must monitor both nitrate and nitrite. An excellent solution to this problem is the enzymatic conversion of nitrate to nitrite by the enzyme nitrate reductase (NaR), followed by quantization of nitrite using Griess Reagent. This kit allows for efficient high-throughput screening of NOS activity in resting cells or cell lysates as well as biological fluids and tissue homogenates. The kit is also ideal for in vitro NOS assays using recombinant purified NOS. All materials necessary to perform the entire assay in a 96-well microplate format are provided with the kit.

For further information about this kit, its application or the procedures in this insert, please contact the Technical Service Team at Eagle Biosciences, Inc. at info@eaglebio.com or at 866-411-8023.

PRINCIPLES OF THE PROCEDURE

NADPH and L-arginine are required for the continual operation of NOS and production of nitric oxide (NO). In aqueous solution, NO rapidly degrades to nitrate and nitrite. Spectrophotometric quantization of nitrite using Griess Reagent is straightforward but does not measure nitrate. This kit employs recombinant nitrate reductase (NaR) for conversion of nitrate to nitrite prior to quantization of nitrite using Griess reagent — thus providing for accurate determination of total NOS activity.



This Nitric Oxide Synthase Ultrasensitive Colorimetric Assay kit can be used to accurately measure as little as 1 pmol/ μ L (~1 μ M) NO produced in aqueous solutions. Very little sample is required (5 to 100 μ L depending on the [NO] in the sample). The completed reaction is read at 540 nm.



MATERIALS PROVIDED

1. **Microtiter Plates** – 1, 96-well low binding, flat-bottom
2. **Reaction Buffer** – 1 bottle containing 50mL of 20mM HEPES 0.5mM EDTA
3. **NADPH Part A** – 1 vial containing 1.0mL of NADP⁺, Glucose 6-Phosphate, L-Arginine **STORE AT -20C**
4. **NADPH Part B** – 1 vial containing 1.0mL of Glucose 6-Phosphate dehydrogenase **STORE AT -20C**
5. **Nitrate Reductase** – 1 vial containing 1U of Nitrate Reductase **STORE AT -20C**
6. **Nitrate Reductase Buffer** – 1.2mL of buffer
7. **Color Reagent #1** – 1 bottle containing 7.0 mL of Sulfanilamide in 3N HCl
8. **Color Reagent #2** – 1 bottle containing 7.0 mL of N-(1-Naphthyl) ethylenediamine dihydrochloride
9. **Nitrite Standard** – 1 vial containing 1.5mL of 500pmol / μ L NaNO₂

MATERIALS NEEDED BUT NOT PROVIDED

1. Reagent grade water: Distilled and deionized.
2. Microplate reader with 540nm filter. (Note: The wavelength of the filter can be 530 to 560 nm, but 540 nm is the absorbance maximum).
3. Precision pipettes ranging from 5 μ L - 1.0 mL and disposable tips. NOTE: If all 96 wells are to be used at one time it is suggested that a multi-channel pipettor be used.
4. Test tubes to dilute the standards and unknowns.
5. Microcentrifuge tubes or microtiter plate for incubation of cell extracts or purified NOS.
6. Centrifuge for microcentrifuge tubes or microtiter plate if incubating in microtiter plate.

STORAGE CONDITIONS

The Ultrasensitive Colorimetric NOS Assay kit is shipped on wet ice and stored at 4-8C. ***Upon arrival, the Nitrate Reductase, NADPH Part A and NADPH Part B should be stored at -20°C until the time of use.*** After initial use of the kit, store each component according to the Materials Provided

PROCEDURAL NOTES

1. Assay can be used with cell lysates or purified NOS for in vitro assays. If using purified NOS it is important the required cofactors are added to the reaction. OBR sells a NOS Cofactor Mix specifically designed for use with this kit under product number **NS70**.
2. If using cell lysates for the assay, it is recommended that they are suspended in phosphate buffered saline (PBS) and protein concentration determined.
3. It is possible for the entire assay to be run in a 96-well microplate if the researcher has a centrifuge equipped with a microplate rotor. The initial incubation should be carried out in the 96-well V-bottom plate while the Color Reaction should be done in the 96-well flat bottom plate.
4. For ease in setting up multiple samples, the Reaction Buffer can be combined with NADPH Part A and Part B and then aliquoted equally to each sample.
5. If the NO concentration in your sample is low, you can increase the sample volume in the Color Reaction to 100 μ L while decreasing the buffer volume.



6. For best results, complete the reading of the plate within 20 minutes.

PREPARATION

1. Nitrate Reductase: Reconstitute with 1ml of Nitrate Reductase Buffer. Vortex briefly every five minutes for fifteen minutes total.
2. Nitrate Standard: Nitrite Standard must first be diluted from 500 μM to 100 μM by adding 240 μL of the Nitrite Standard to 960 μL of Reaction Buffer. This should be labeled Diluted Standard and stored on ice until used. The standard curve is then created by further diluting the Diluted Standard according to Table 2 below.

Table 2: Preparation of the Standard Curve.

Standard	Final Concentration (μM)	ddH ₂ O (μL)	Diluted Standard (μL)
B ₀	0	1000	0
S ₁	0.5	995	5
S ₂	1	990	10
S ₃	5	950	50
S ₄	10	900	100
S ₅	25	750	250
S ₆	50	500	500
S ₇	100	0	250

NOTE: Standard solutions may be stored at 4°C for later use.

PROCEDURE

In V-well microplate or microcentrifuge tube

1. Add 40–500 μg of protein from lysates or 0.2 - 1.0 Unit of recombinant or purified NOS in a volume of 30 μL to a tube or well.
2. Add 200 μL Reaction Buffer
3. Add 10 μL of NADPH Part A
4. Add 10 μL of NADPH Part B
5. If using purified NOS, add 10 μL of the NS70 – NOS Cofactor Mix
6. Mix and incubate for 1 - 6 hours at 37°C.
7. Chill on ice for 5 minutes.
8. Add 10 μL of the reconstituted Nitrate Reductase to each sample, vortex tube or tap plate to mix, and incubate for 20 minutes at room temperature.
9. Centrifuge at 12,500 rpm for 5 minutes at 4°C.

In flat-bottom microtiter plate

10. Add 100 μL of Standards or Samples to the microplate in duplicate. Depending on [NO], the samples may need to be diluted in Reaction Buffer. See Scheme 1 for a sample plate layout.
11. Add 50 μL Color Reagent #1 and shake briefly.



12. Add 50 μL Color Reagent #2. Shake for 5 minutes at room temperature.
13. Read absorbance values at 540 nm in Microtiter plate reader.

Table 3: Sample microplate layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	B ₀	B ₀	U ₁	U ₁	U ₉	U ₉	U ₁₇	U ₁₇	U ₂₅	U ₂₅	U ₃₃	U ₃₃
B	S ₁	S ₁	U ₂	U ₂	U ₁₀	U ₁₀	U ₁₈	U ₁₈	U ₂₆	U ₂₆	U ₃₄	U ₃₄
C	S ₂	S ₂	U ₃	U ₃	U ₁₁	U ₁₁	U ₁₉	U ₁₉	U ₂₇	U ₂₇	U ₃₅	U ₃₅
D	S ₃	S ₃	U ₄	U ₄	U ₁₂	U ₁₂	U ₂₀	U ₂₀	U ₂₈	U ₂₈	U ₃₆	U ₃₆
E	S ₄	S ₄	U ₅	U ₅	U ₁₃	U ₁₃	U ₂₁	U ₂₁	U ₂₉	U ₂₉	U ₃₇	U ₃₇
F	S ₅	S ₅	U ₆	U ₆	U ₁₄	U ₁₄	U ₂₂	U ₂₂	U ₃₀	U ₃₀	U ₃₈	U ₃₈
G	S ₆	S ₆	U ₇	U ₇	U ₁₅	U ₁₅	U ₂₃	U ₂₃	U ₃₁	U ₃₁	U ₃₉	U ₃₉
H	S ₇	S ₇	U ₈	U ₈	U ₁₆	U ₁₆	U ₂₄	U ₂₄	U ₃₂	U ₃₂	U ₄₀	U ₄₀



S₀ = Blanks: Add buffer in place of standard or sample.

U₁₋₄₀ = Unknown sample to be assayed.

CALCULATIONS

1. Subtract the average O.D. value of the blank wells (S₀) from all other pairs of wells.
2. Average the O.D. values for each pair of duplicate wells.
3. Plot a standard curve using the average O.D. value for each standard value versus the concentration of standard. If samples are less than 10 μM , plot standards B₀, S₁, S₂, S₃, and S₄. For samples over 10 μM , plot standards S₄, S₅, S₆, and S₇.
4. Determine the concentration of each unknown by interpolation from the standard curve.
5. Samples can be compared by determining μmoles of NO produced / μg protein / unit time.

Figure 1: Typical Standard Curve





WARNINGS AND PRECAUTIONS

1. **DO NOT** use components beyond the specified expiration date.
2. Universal precautions should be employed with handling of the components provided in this kit.
3. Use aseptic techniques when opening and removing reagents from vials and bottles.

REFERENCES

Ghigo, D., et al. (2006). *Nitric Oxide* **15**:148-153
Schmidt, H.H., et al. (1995). *Biochemica* **2**:22-23

WARRANTY INFORMATION

Eagle Biosciences, Inc. warrants its Product(s) to operate or perform substantially in conformance with its specifications, as set forth in the accompanying package insert. This warranty is expressly limited to the refund of the price of any defective Product or the replacement of any defective Product with new Product. This warranty applies only when the Buyer gives written notice to the Eagle Biosciences within the expiration period of the Product(s) by the Buyer. In addition, Eagle Biosciences has no obligation to replace Product(s) as result of a) Buyer negligence, fault, or misuse, b) improper use, c) improper storage and handling, d) intentional damage, or e) event of force majeure, acts of God, or accident.

Eagle Biosciences makes no warranties, either expressed or implied, except as provided herein, including without limitation thereof, warranties as to marketability, merchantability, fitness for a particular purpose or use, or non-infringement of any intellectual property rights. In no event shall the company be liable for any indirect, incidental, or consequential damages of any nature, or losses or expenses resulting from any defective product or the use of any product. Product(s) may not be resold, modified, or altered for resale without prior written approval from Eagle Biosciences, Inc.

For further information about this kit, its application, or the procedures in this kit, please contact the Technical Service Team at Eagle Biosciences, Inc. at info@eaglebio.com or at 866-411-8023.