## **GroPep Bioreagents**



# **Barramundi IGF-I antiserum (Rabbit)**

#### Description

Barramundi IGF-I antiserum was prepared by immunising rabbits with recombinant barramundi IGF-I. It has 100% cross-reactivity with barramundi IGF-I and other species of fish when determined by immunoassay. It can be used as a ubiquitous antiserum for the determination of plasma or serum IGF-I in a range of fish species.

**References:** Degger B et al (2000) General and Comparative Endocrinology 117, 395-

403

**Immunogen:** Barramundi (*Lates calcarifer*) IGF-I

**Source:** Rabbit antiserum

**Specificity:** Cross-reacts 100% with Barramundi IGF-I, Bream IGF-I (*Pagarus* 

auratus), Salmon and Trout IGF-I (Onchorynchus and Salmo species) Tilapia IGF-I ( $Oreochromis\ mossambicus$ ) and Tuna IGF-I ( $Thunnus\ maccoyii$ ). Cross-reactivity to human IGF-I and IGF-II is < 1%. No

cross reactivity to salmon insulin.

**Appearance:** Lyophilized powder

**Storage/Stability:** At least 2 years at 2 - 4°C (lyophilized).

After reconstitution store at -20°C or -80°C

Avoid freeze-thaw cycles

Reconstitution: Dissolve the vial contents in 500µl RIA buffer (30mM sodium

phosphate pH 7.5 containing 10mM EDTA, 0.2g/l protamine sulphate, 0.5ml Tween 20, 0.25g/l sodium azide) to give an

antiserum stock of 1:50.

**Application and Titre:** Each application and titre should be determined in house but as a

guideline - Determination of IGF-I or IGF-II in a Range of Species by Radioimmunoassay (RIA). At a dilution of 1:30,000 the antiserum binds 20% of freshly prepared [ $^{125}$ I]- barramundi IGF-I tracer ( $^{80}$ Ci/g). Half-maximal inhibition (ED50) of binding is produced by  $^{80}$ 

1.7 ng/ml of barramundi IGF-I.

#### **Product Codes and Pricing**

**Barramundi IGF-I** 

antiserum 10 ul PAF1

Related Products: Barramundi IGF-I

Bream IGF-I

Salmon/Trout IGF-I

Tuna IGF-I

### \*\*NOT FOR USE IN HUMANS\*\*

**GroPep Bioreagents Pty Ltd** 51 West Thebarton Road Thebarton SA 5031

Australia

Postal Address: PO Box 10065

Adelaide Business Centre

SA 5000 Australia Email:

info@gropep.com

Internet: www.gropep.com

**ABN** 93 147 032 166

Telephone: +61 8 7222 1051