



# APIS DPYD qPCR Kit

The APIS DPYD qPCR Kit is designed for fast and efficient detection of seven *DPYD* alleles. Each reaction identifies both wildtype (WT) and variant alleles, directly from blood or extracted genomic DNA (gDNA). The kit's rapid turnaround time, cost-effectiveness, and easy integration into existing qPCR workflows makes it a reliable and accessible option for a variety of research applications.

## Advantages of the APIS DPYD qPCR Kit



Wide target coverage – detects seven DPYD variants



Direct workflow – compatible with both 1% EDTA Whole Blood and gDNA



Reliable HapB3 detection – targets key functional HapB3 variant



Accessible – compatible with most qPCR instruments

### **Analytical Sensitivity**

The analytical sensitivity of the kit was evaluated by testing homozygous (wildtype) and heterozygous samples and defined as the minimum amount of DNA necessary to correctly genotype 95% of the samples. The APIS DPYD qPCR Kit can detect seven *DPYD* variants in as little as 0.11ng gDNA.

Table 1: DPYD variants and analytical sensitivity of the APIS DPYD qPCR Kit

| Variant (NM_000110.4) | rsID        | Sensitivity (ng per reaction) |
|-----------------------|-------------|-------------------------------|
| c. 1905+1G>A          | rs3918290   | 0.11                          |
| c.1679T>G             | rs55886062  | 0.04                          |
| c.1129-5923C>G        | rs75017182  | 0.11                          |
| c.557A>G              | rs115232898 | 0.04                          |
| c. 2846A>T            | rs67376798  | 0.04                          |
| c.868A>G              | rs146356975 | 0.04                          |
| c.2279C>T             | rs112766203 | 0.04                          |

# APIS DPYD qPCR Kit shows excellent concordance with alternative research genotyping methods

gDNA extracted from 29 specimens with known *DPYD* variant status (characterised via NGS, microarrays or qPCR) and one synthetic DNA sample were assessed with the APIS DPYD qPCR Kit.

Table 2: Concordance of the APIS DPYD qPCR Kit with other methods

| Variant<br>(NM_000110.4) | N° Homozygous WT<br>specimens | N° Heterozygous<br>specimens | N° Homozygous<br>Variant specimens | N° of correctly<br>genotyped specimens |
|--------------------------|-------------------------------|------------------------------|------------------------------------|--|
| c.1905+1G>A              | 24                            | 4                            | 2                                  | 30/30                                  |
| c.1679T>G                | 27                            | 2                            | 1                                  | 30/30                                  |
| c.1129-5923C>G           | 23                            | 5                            | 2                                  | 30/30                                  |
| c.557A>G                 | 24                            | 5                            | 1                                  | 30/30                                  |
| c.2846A>T                | 24                            | 5                            | 1                                  | 30/30                                  |
| c.868A>G                 | 26                            | 3                            | 1                                  | 30/30                                  |
| c.2279C>T                | 24                            | 5                            | 1                                  | 30/30                                  |

#### APIS DPYD qPCR Kit Specifications

The kit has undergone testing as part of APIS unique RUO assay development pipeline, including both *in silico* and wet laboratory testing. As a result, researchers can depend on the kit for consistent and high-quality detection of *DPYD* variants. A positive and negative control (PC and NTC), which monitor assay set up and assay performance, are also supplied with the kit.

| Format               | 11 Tubes: 7 PP mixes (Primer-Probes), 2 Enzyme mixes, 1 Positive Control, 1 Negative Control |  |  |
|----------------------|--|--|--|
| PCR Method           | qPCR   |  |  |
| Turnaround Time      | < 2 hours per PCR plate  |  |  |
| Hands-on Time        | < 1 hours per PCR plate  |  |  |
| Intended Sample Type | Extracted gDNA (e.g. from whole blood) or 1% EDTA whole blood                                |  |  |
| Instrument           | qPCR instrument  |  |  |
| Shelf Life           | 12 months stored at -25°C to -15°C   |  |  |

# APIS DPYD qPCR Kit Ordering Information

| Assay Name         | Test Type | Kit Size                 | Catalogue Number | Price                  |
|--------------------|-----------|--------------------------|------------------|------------------------|
| APIS DPYD qPCR Kit | RUO       | 24 samples plus controls | 03501            | Available upon request |

All scientific and analytical information in this brochure comes from the APIS DPYD qPCR Kit Handbook ART0180.

The APIS DPYD qPCR Kit is intended for Research Use Only (RUO). Not for use in diagnostic procedures.