

HLA-A*24:07 easYmer kit

Catalog#	1023-01
Size	20, 50, 150, or 500 tests
easYmer	HLA-A*24:07 - peptide receptive, biotinylated in Tris/Maleate pH 7 with 30% Glycerol
Folding buffer	Tris/Maleate pH 7
Peptide	SYMSTFPLF, a stable HLA-A*24:07 binder Positive control for evaluation analysis of peptide-HLA folding
Application	easYmers can be used to generate monomers with your choice of peptide. The monomers can easily be tetramerized with fluorophore conjugated streptavidin and used to analyse T cells by flow cytometry. The easYmer reagent can also be used to evaluate specific peptide-HLA I binding.
	For Research Use Only (RUO)