New Tools for Diagnostics in Autoimmunity



Simply innovative diagnostics



Fully automated Dot/LINE Processor



Technical features:

- Kit processing and evaluation in one with minimal user intervention
- Simultaneous performance of up to 24 different tests
- Barcode identification of test strips and reagent cartridges
- Digitalized results
- LIMS ready, results on paper optional
- Walk-away-system, just load and press start
- Very low maintenance, no hydraulic pumps, multiple washings unnecessary
- Small footprint





DotDiver2.0 assays available for:

- Celiac disease / Crohn's disease
- Pernicious anemia
- Neuropathies
- Anti-phospholipid syndrome
- Systemic lupus erythematosus
- Vasculitis
- Nephritis
- Myositis
- Scleroderma
- Autoimmune hepatitis

Dot/LINE Immunoassays: Automated processing and interpretation


Test procedure



Insertion of strips and reagents



Incubation and washing



Strip drying and read out

Technical Parameters



Results (example)

Dimensions Weight	34 x 42 x 36.5 cm (L x W x H) 15 kg
Assay combinations per run	Maximum 24
Sample capacity	1 - 24
Method	Antigen coated strips ready to use reagent cartridge
Assay time	1.5 h
Software	Integrated evaluation software
Safety	Barcode identification of strips and cartridges
Results	LIMS or paper print out (optional)
Indications	Systemic inflammatory rheumatic diseases, gastroenterology, neuropathies, anti phospholipid syndrome, vasculitis, nephritis, myositis, scleroderma, autoimmune hepatitis
	More than 20 different test kits available
Order-No.	5075





GA Generic Assays GmbH Ludwig-Erhard-Ring 3 15827 Dahlewitz/Berlin Germany Phone: +49 33708 9286 0 Fax: +49 33708 9286 50 info@genericassays.com www.genericassays.com