

Section 1: Product and company information

Product identifier

Product name	Chitozen 5-coverslip Standard Kit
Substance name	N/A
Reach No.	A registration number is not available for this substance as the substance or its uses are exempted from registration, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline

HTS Code 3822.00

Relevant identified uses of the substances or mixture and uses advised against

Relevant Identified uses

Main use category	Laboratory use
Professional use spec	R&D use only
Identified used	Laboratory chemicals

Details of the supplier of the safety data sheet

Company	Idylle 68 boulevard de Port-Royal 75005 Paris France
Technical phone	+33184172227
Email address	contact@idylle-labs.com

Emergency phone number

France	112
Other countries	ORFILA (INRS) : + 33 (0)1 45 42 59 59

Section 2: Technical information

Kit and product description

- Each kit is composed of **5 chitosan-coated coverslips** and packaged in a shockproof box in which up to 5 standard units can be stored.
- Each Chitozen Standard coverslip is composed of the **chitosan-coated coverslip**, standing against a **protective slide** on which the product label stands. The chitosan-coated coverslip is originally linked to the protective slide by two removable adhesives. These two removable adhesives are to be peeled off prior to use (see protocol). The chitosan-coated side of the coverslip is the one facing up (see scheme below).
- Dimensions of the chitosan-coated coverslip (w x l x h): 75 x 25 x 0.175 mm³



Shipping, storage and expiration

- Each kit of **5 chitosan-coated coverslips** is shipped at room temperature and shielded from direct sunlight light in the shockproof box.
- The **chitosan-coated coverslips** should be stored at room temperature, away from moisture and shielded from direct sunlight, in the shockproof box they are originally packed.
- The **chitosan-coated coverslips** can be stored and used up to 12 months after production. Production year and month is indicated by the batch code

According to Regulation (EC) No. 1907/2006(REACH) with its amendment Regulation (EC) No. 453/2010-830/2015

YYMM (where YY stands for the last two digits of the year and MM for the month).

distributed in the US/Canada by:

EAGLE BIOSCIENCES, INC.

20A NW Blvd, Suite 112 Nashua, NH 03063

Phone: 617-419-2019 FAX: 617-419-1110

www.EagleBio.com • info@eaglebio.com



EAGLE
BIOSCIENCES