



EAGLE
BIOSCIENCES

Anti-Protein Kinase C eta (η)

Catalog Number: KCE00-A100

Aliquot: 100 μ L

Storage: -20°C

PRODUCT SPECIFICATION SHEET

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C- η (PKC η).
MOLECULAR WEIGHT:	67.7 kDa Calculated
CONCENTRATION:	Not Determined
PRESENTATION:	Antiserum with 0.05% Sodium Azide
STORAGE:	-20°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
IMMUNOGEN:	KLH conjugated synthetic peptide corresponding to amino acids 666- 683 of PKC η . Sequence: ⁶⁶⁶ N-Q-D-E-F-R-N-F-S-Y-V-S-P-E-L-Q-P ⁶⁸³
PURIFICATION:	Not Purified
APPLICATIONS:	Western Blotting: 1:1,000 dilution using ECL detection technology.
SPECIES SPECIFICITY:	Human, mouse and rat among other mammalian species.

Product Developed and Manufactured in the USA

EAGLE BIOSCIENCES, INC.
20A NW Blvd., Suite 112
Nashua, NH 03063

Phone: 617-419-2019
Fax: 617-419-1110
Email: info@EagleBio.com
www.EagleBio.com