



EAGLE
BIOSCIENCES

Anti-Protein Kinase C β I

Catalog Number: PK100-A100

Aliquot: 100 μ L

Storage: -70°C

PRODUCT SPECIFICATION SHEET

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C β I (PKC β I).
MOLECULAR WEIGHT:	76.9 kDa Calculated
CONCENTRATION:	Not Determined
ISOTYPE:	NA
STORAGE BUFFER:	Phosphate Buffered Saline; 0.05% Sodium Azide
STORAGE:	-70°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
IMMUNOGEN:	KLH conjugated synthetic peptide corresponding to amino acids 661- 671 of PKC β I. Sequence: ⁶⁶¹ S-Y-T-N-P-E-F-V-I-N-V ⁶⁷¹
PURIFICATION:	NA
APPLICATIONS:	Western Blotting: 1:1000 Dilution
SPECIES SPECIFICITY:	Human, mouse and rat among other mammalian species.
REFERENCE:	JBC (1993) 268 :20110-20115

Product Developed and Manufactured in the USA

EAGLE BIOSCIENCES, INC.
20A NW Blvd., Suite 112
Nashua, NH 03063

Phone: 617-419-2019
Fax: 617-419-1110
Email: info@EagleBio.com
www.EagleBio.com