

Product data



(IVTDFSVIK)/HLA-A*11:01

Catalog#	1018-03
Size	50, 150, or 500 tests
Organism	EBV
Epitope	EBNA-4 416-424 (IEDB#: 29466)
HLA	Recombinant HLA-A*11:01 with a C-terminal biotinylation.
Formulations	peptide-HLA-I monomer, biotinylated peptide-HLA-I tetramer with PE or APC label
Buffer	Tris/Maleate PBS pH 7 containing 5% Glycerol
Application	The tetramers can be used in flow cytometric analysis to characterise and quantify epitope specific CD8+T cells. The biotinylated monomer can be made into tetramers with optional labels.

For Research Use Only (RUO)