

Product data



(IYGIFQSTF)/HLA-A*24:02

Catalog#	1020-04
Size	50, 150, or 500 tests
Organism	Yellow Fever virus
Epitope	NS3 24-32 (IEDB#: 232000)
HLA	Recombinant HLA-A*24:02 with a C-terminal biotinylation.
Formulations	peptide-HLA-I monomer, biotinylated peptide-HLA-I tetramer with PE or APC label
Buffer	Tris/Maleate pH 7 containing 5% Glycerol
Application	The tetramers can be used in flow cytometric analysis to characterise and quantify epitope specific CD8+T cells. The biotinylated monomer can be made into tetramers with optional labels.

For Research Use Only (RUO)