

(RPPIFIRRL)/HLA-B*07:02

Catalog#	1048-02
Size	50, 150, or 500 tests
Organism	EBV
Epitope	EBNA-3A 247-255 (IEDB#: 55252)
HLA	Recombinant HLA-B*07:02 with a C-terminal biotinylation.
Formulations	peptide-HLA-I monomer, biotinylated peptide-HLA-I tetramer with PE or APC label
Buffer	Tris/Maleate pH 7 containing 5% Glycerol
Application	The tetramers can be used in flow cytometric analysis to characterise and quantify epitope specific CD8+T cells. The biotinylated monomer can be made into tetramers with optional labels.

For Research Use Only (RUO)