

### (TPRVTGGGAM)/HLA-B\*07:02

<b>Catalog#</b>	1048-03
<b>Size</b>	50, 150, or 500 tests
<b>Organism</b>	HCMV
<b>Epitope</b>	pp65 417-426 (IEDB#: 65748)
<b>HLA</b>	Recombinant HLA-B*07:02 with a C-terminal biotinylation.
<b>Formulations</b>	peptide-HLA-I monomer, biotinylated peptide-HLA-I tetramer with PE or APC label
<b>Buffer</b>	Tris/Maleate pH 7 containing 5% Glycerol
<b>Application</b>	The tetramers can be used in flow cytometric analysis to characterise and quantify epitope specific CD8+T cells. The biotinylated monomer can be made into tetramers with optional labels.

**For Research Use Only (RUO)**