

Product data



(FPKTTNGCSQA)/HLA-B*55:01

Catalog#	1103-02
Size	50, 150, or 500 tests
Organism	HCMV
Epitope	IE1 221-231 (IEDB#: 17321)
HLA	Recombinant HLA-B*55:01 with a C-terminal biotinylation.
Formulations	peptide-HLA-I monomer, biotinylated peptide-HLA-I tetramer with PE or APC label
Buffer	Tris/Maleate pH 7 containing 5% Glycerol
Application	The tetramers can be used in flow cytometric analysis to characterise and quantify epitope specific CD8+T cells. The biotinylated monomer can be made into tetramers with optional labels.

For Research Use Only (RUO)