

Product data



(KGIQIYTRNHEVKS)/HLA-DRA1*01:01/HLA-DRB1*07:01

Catalog#	2009-02
Size	50, 150, or 500 tests
Organism	HCMV
Epitope	IE2 408-422 (IEDB#: 230003)
HLA	Recombinant HLA-DRA1*01:01 / recombinant HLA-DRB1*07:01 with a C-terminal biotinylation.
Formulations	peptide-HLA-II monomer, biotinylated peptide-HLA-II tetramer with PE or APC label
Buffer	PBS pH 7.4 containing 5% Glycerol
Application	The tetramers can be used in flow cytometric analysis to characterise and quantify epitope specific CD4+T cells. The biotinylated monomer can be made into tetramers with optional labels.

For Research Use Only (RUO)