

## Pre-immune normal rabbit IgG

### Description

This immunoglobulin G (IgG) fraction was purified from normal rabbit serum obtained from animals with no history of immunisation. It can be used as a negative control for GroPep Bioreagents' range of purified rabbit IgG antibodies.

**Source:** Affinity purified from normal rabbit serum

**Appearance:** Lyophilized powder

**Storage/Stability:** At least 2 years at 2 - 4°C (lyophilized).  
After reconstitution store at -20°C or -80°C  
Avoid freeze-thaw cycles

**Reconstitution:** Dissolve the vial contents in 1 ml sterile water and then dilute in a buffer suitable for the application.

### Product Codes and Pricing

**Normal Rabbit IgG**      **100 µg**      **NRI01**

**Related Products:** Pre-immune normal rabbit serum

**\*\*NOT FOR USE IN HUMANS\*\***

**GroPep Bioreagents Pty Ltd**  
51 West Thebarton Road  
Thebarton SA 5031  
Australia

**ABN** 93 147 032 166

**Telephone:** +61 8 7222 1051

**Postal Address:**  
PO Box 10065  
Adelaide Business Centre  
SA 5000  
Australia

**Email:**  
[info@gropep.com](mailto:info@gropep.com)

**Internet:**  
[www.gropep.com](http://www.gropep.com)