

Pre-immune normal rabbit serum

Description

This normal rabbit serum was obtained from animals with no history of immunisation. It can be used as a negative control for GroPep Bioreagents' range of antisera raised in rabbits.

Source: Normal rabbit serum

Appearance: Lyophilized powder

Storage/Stability: At least 2 years at 2 - 4°C (lyophilized).
After reconstitution store at -20°C or -80°C
Avoid freeze-thaw cycles

Reconstitution: Dissolve the vial contents in 200 µl sterile water and then dilute in a buffer suitable for the application.

Product Codes and Pricing

Normal rabbit serum **200 µl** **NRS01**

Related Products: Pre-immune normal rabbit IgG

****NOT FOR USE IN HUMANS****

GroPep Bioreagents Pty Ltd
51 West Thebarton Road
Thebarton SA 5031
Australia

ABN 93 147 032 166

Telephone: +61 8 7222 1051

Postal Address:
PO Box 10065
Adelaide Business Centre
SA 5000
Australia

Email:
info@gropep.com

Internet:
www.gropep.com