GroPep Bioreagents



Salmon/Trout IGF-II Polyclonal Antibody (Rabbit)

Description

Salmon/Trout IGF-II polyclonal antibody was prepared by immunising rabbits with a synthetic peptide epitope corresponding to amino acid residues 1-9 of the N-terminus of mature Atlantic salmon (Salmo salar) IGF-II. The antiserum was harvested and the IgG fraction partially purified by hydrophobic charge induction chromatography. The antibody is highly specific for IGF-II from Atlantic salmon and Rainbow trout (Oncorhynchus mykiss) in an enzyme-linked immunoassay. It has a cross-reactivity of ~ 5% to salmon insulin and < 0.01% to salmon IGF-I and salmon GH.

Reference: Wilkinson R et al (2004) Comparative Biochemistry and Physiology B 139,

193 - 201.

Immunogen: Derived from UniProtKB/Swiss-Prot sequence ID A7LBA7 Atlantic

salmon IGF-II. The peptide sequence is amino-acids 11-19 (IGF-II₁₋₉)

conjugated to diphtheria toxoid.

Source: HCIC purified IgG from rabbit antiserum (100 µl). Each vial contains

approximately 520 µg immunoglobulin.

Specificity: 5.6 % cross-reactivity with salmon insulin

> < 0.1 % cross-reactivity with salmon/trout IGF-I (Oncorhynchus kisutch & Oncorhynchus mykiss), salmon/trout growth hormone, tuna IGF-I (Thunnus maccoyii), tuna IGF-II, bream IGF-I (Pagrus auratus), bream growth hormone (Acanthopagrus butcheri), barramundi IGF-I (Lates calcarifer), human IGF-I, human IGF-II, rat

IGF-I, rat IGF-II, rat growth hormone.

Appearance: Lyophilized powder

Storage/Stability: At least 2 years at 2 - 4°C (lyophilized).

After reconstitution store at -20°C or -80°C

Avoid freeze-thaw cycles

Reconstitution: Dissolve the vial contents in 500µl buffer suitable for enzyme-linked

immunoassay applications.

Application and Titre: Each application and titre should be determined in house but as a

guideline:

Enzyme immunoassay 1:1,000 Immunoblotting 1:1,000

Product Codes and Pricing

Salmon/Trout IGF-II Polyclonal antibody 100 µl PAAA1

Related Products: Salmon/Trout IGF-II

NOT FOR USE IN HUMANS

GroPep Bioreagents Pty Ltd 51 West Thebarton Road Thebarton SA 5031

Australia

Postal Address: PO Box 10065

Adelaide Business Centre

SA 5000 Australia info@gropep.com

Internet: www.gropep.com

ABN 93 147 032 166

Telephone: +61 8 7222 1051