



GroPep Bioreagents

Rabbit IgG antiserum (Goat)

Description

This anti-rabbit immunoglobulin G (IgG) antiserum is harvested from goats which were immunised with affinity purified rabbit IgG. The product can be used as a secondary antibody in a number of immunoassay protocols which employ a rabbit IgG as the primary antibody.

Immunogen:	Affinity purified rabbit IgG
Source:	Goat antiserum
Appearance:	Lyophilized powder
Storage/Stability:	At least 2 years at 2 - 4°C (lyophilized). After reconstitution store at -20°C or -80°C Avoid freeze-thaw cycles
Reconstitution:	Dissolve in a buffer suitable for the application.

Product Codes and Pricing

Rabbit IgG antiserum	625 µl	PSA1
-----------------------------	---------------	-------------

Related Products:	Pre-immune normal rabbit serum Pre-immune normal rabbit IgG
--------------------------	--

****NOT FOR USE IN HUMANS****

GroPep Bioreagents Pty Ltd
51 West Thebarton Road
Thebarton SA 5031
Australia

ABN 93 147 032 166

Telephone: +61 8 7222 1051

Postal Address:
PO Box 10065
Adelaide Business Centre
SA 5000
Australia

Email:
info@gropep.com

Internet:
www.gropep.com