

Rat IGF-I

Description

Rat insulin-like growth factor-I (IGF-I) is a single chain 70 amino acid polypeptide that shares structural homology with human and rat IGF-II. Only three amino acids are different compared to the sequence of human IGF-I. Rat IGF-I has been shown to stimulate growth in L6 myoblasts and other cell lines. Receptor Grade rat IGF-I is a high quality research reagent particularly for use in studies on cellular and animal growth, for use in specific assays for rat IGF-I and other general research purposes.

References

Shimatzu A and Rotwein P. (1987) Journal of Biological Chemistry **262**, 7894-7900
Upton Z *et al* (1996) Journal of Endocrinology **149**, 379-387

Source:	Produced in <i>E.coli</i> .
Molecular Weight:	7687 daltons
Purity:	>95 % (by HPLC analysis)
N-terminal sequence:	Analysis of 5 residues
Biological Activity:	IGF type I receptor binding assay (ED ₅₀ < 10 ng/ml) Stimulation of protein synthesis in rat L6 myoblasts (ED ₅₀ < 30 ng/ml)
Endotoxin:	< 0.1 EU/μg
Appearance:	White powder freeze-dried from 0.1M acetic acid and stored under nitrogen at a slight vacuum.
Storage/Stability:	At least 2 years at 2-4°C (as a freeze dried product).
Reconstitution:	<i>Handling of GroPep IGF-I, IGF-II and IGF analogues</i>
Detection:	<i>Procedure for the Iodination of IGF peptides</i> <i>Determination of IGF-I or IGF-II in a range of species by Radioimmunoassay (RIA)</i>

Product Codes:

Rat IGF-I	100 μg	RU100
	1 mg	RM001

Related Products:

Rat IGF-II
Rat IGF-I antiserum (Rabbit)
IGF-I antibody for immunohistochemistry
Rat IGF-I monoclonal antibody

****NOT FOR USE IN HUMANS****

GroPep Bioreagents Pty Ltd
51 West Thebarton Road
Thebarton SA 5031
Australia

ABN 93 147 032 166

Telephone: +61 8 7222 1051

Postal Address:
PO Box 10065
Adelaide Business Centre
SA 5000
Australia

Email:
info@gropep.com

Internet:
www.gropep.com