

GroPep Bioreagents

Human IGF type 1 receptor 26-3 monoclonal antibody (alpha chain)

Description

Human IGF type 1 receptor antibody 26-3 is a mouse monoclonal antibody (IgG1) which recognizes an epitope within amino acids 283-440 or 586-908 of the extracellular domain (α -chain) of the human IGF type 1 receptor. The antibody was prepared by immunising mice with a preparation of 3T3 cells which over-express the IGF type 1 receptor. The antibody cross reacts with the human IGF type 1 receptor but not the human or rat insulin receptor.

References: Soos MA *et al* (1992) Journal of Biological Chemistry **267**, 12955-12963

Immunogen: 3T3 cells over-expressing the IGF type 1 receptor

Source: Produced in hybridoma cells derived from the NS-1 myeloma cell line. Purified from conditioned medium by Protein A affinity chromatography

Specificity: Cross-reacts with the human IGF type 1 receptor but not the human or rat insulin receptor

Appearance: Lyophilized white powder

Storage/Stability: At least 2 years at 2 - 4°C (lyophilized).
After reconstitution store at -20°C or -80°C
Avoid freeze-thaw cycles

Reconstitution: Dissolve in 200 μ l phosphate buffered saline pH 7.4

Application and Titre: Immunoprecipitation (1:5,000)
Enzyme immunoassay (1:5,000)

Product Codes and Pricing

Human IGF type I receptor monoclonal antibody 26-3	200 μg	MAF1
---	------------------------------	-------------

Related Products: Human IGF type 1 receptor monoclonal antibody CT-1;
Human IGF type 1 receptor monoclonal antibody 24-55;
Human IGF type 1 receptor monoclonal antibody 16-13

****NOT FOR USE IN HUMANS****

GroPep Bioreagents Pty Ltd
51 West Thebarton Road
Thebarton SA 5031
Australia

ABN 93 147 032 166

Telephone: +61 8 8152 9360

Postal Address:
PO Box 10065
Adelaide Business Centre
SA 5000
Australia

Email:
info@gropep.com

Internet:
www.gropep.com