



**EAGLE**  
BIOSCIENCES

## Anti-6-Keto-Prostaglandin F1 $\alpha$

Catalog Number: 6KP00-A1000

Aliquot: 1000  $\mu$ l

Storage: -80°C

### PRODUCT SPECIFICATION SHEET

<b>DESCRIPTION:</b>	This rabbit polyclonal antibody is specific to 6-Keto-Prostaglandin F1 $\alpha$ and is provided as antiserum.
<b>MOLECULAR WEIGHT:</b>	36 kDa
<b>CONCENTRATION:</b>	Not determined
<b>STORAGE BUFFER:</b>	PBS; 0.1% Sodium Azide
<b>STORAGE:</b>	-80°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	6-Keto-Prostaglandin F1 $\alpha$
<b>SENSITIVITY:</b>	150 fmol/assay (IC <sub>50</sub> )
<b>PURIFICATION:</b>	Not purified
<b>APPLICATIONS:</b>	Suitable for use in RIA. May also have utility in ELISA.
<b>SPECIES SPECIFICITY:</b>	Not species specific. May not exhibit utility when using in a rabbit experimental model.
<b>REFERENCE:</b>	Levine, L. Measurement of Arachidonic Acid Metabolites by Radioimmunoassay. <u>In</u> : <i>Manual of Clinical Laboratory Immunology</i> (3rd Edition). J.L. Fahey, H. Friedman & N.L. Rose (eds). Amer. Soc. Microbiol. pp. 685-691 (1986).

*Product Developed and Manufactured in the USA*

EAGLE BIOSCIENCES, INC.  
20A NW Blvd., Suite 112  
Nashua, NH 03063

Phone: 617-419-2019  
Fax: 617-419-1110  
Email: [info@EagleBio.com](mailto:info@EagleBio.com)  
[www.EagleBio.com](http://www.EagleBio.com)