

DATA SHEET

Anti human C4d Antibody (C4dpAb)

Cat. no.: BI- RC4D

For research use only

Source: rabbit polyclonal

Quantity: 250 µl/vial

Form: liquid IgG fraction, purified by protein G chromatography

Immunogen: synthetic peptide (15-mer peptide corresponding to aa 1242-1256 of C4) specific for human complement split product C4d coupled to KLH (keyhole limpet haemocyanine).

Specificity: The antibody will identify human complement split product C4d

Application: The antibody has been tested for use on paraffin sections and frozen sections of human tissue. Working protocol available.
Recommended dilution for immunohistochemistry: 1:30 – 1:50
Pretreatment of slides by pressure cooking (5 minutes at +125°C in TRIS/EDTA buffer pH 8.5) is highly recommended for antigen retrieval.

Stability: Storage and transport temperature: +4°C.
The antibody is stable at +4°C up to the expiry date.

Shelf life: Expiry date is stated on the label.


Literature:

1. Regele H et al. "Endothelial C4d deposition is associated with inferior kidney allograft outcome independently of cellular rejection". Nephrol Dial Transplant 2001, 16:2058-66

2. Böhmig GA et al. C4d-positive acute humoral renal allograft rejection: effective treatment by immunoadsorption. J Am Soc Nephrol, 2001 Nov;12(11):2482-9

3. Böhmig GA et al. Capillary C4d deposition in kidney allografts: a specific marker of alloantibody-dependent graft injury. J Am Soc Nephrol. 2002 Apr;13(4):1091-9.

4. Regele H et al. „Capillary deposition of complement split product C4d in renal allografts is associated with basement membrane injury in peritubular and glomerular capillaries: a contribution of humoral immunity to chronic allograft rejection". J Am Soc Nephrol. 2002, 13:2371-80




rev.no: 120229 (replacing 100705))

 BIOMEDICA




Biomedica Medizinprodukte GmbH & Co KG, A-1210 Wien, Divischgasse 4
Tel. +43/1/291 07 45, Fax +43/1/291 07 85, E-mail export@bmgrp.com

Microphotographs of Anti human C4d Antibody (C4dpAb)




paraffin embedded tissue sections



Cryosections

Immunofluorescence Microphotographs of Anti human C4d Antibody (C4dpAb)



Dilution 1:20

 BIOMEDICA

Warranty Information

Eagle Biosciences, Inc. warrants its Product(s) to operate or perform substantially in conformance with its specifications, as set forth in the accompanying package insert. This warranty is expressly limited to the refund of the price of any defective Product or the replacement of any defective Product with new Product. This warranty applies only when the Buyer gives written notice to the Eagle Biosciences within the expiration period of the Product(s) by the Buyer. In addition, Eagle Biosciences has no obligation to replace Product(s) as result of a) Buyer negligence, fault, or misuse, b) improper use, c) improper storage and handling, d) intentional damage, or e) event of force majeure, acts of God, or accident.

Eagle Biosciences makes no warranties, either expressed or implied, except as provided herein, including without limitation thereof, warranties as to marketability, merchantability, fitness for a particular purpose or use, or non-infringement of any intellectual property rights. In no event shall the company be liable for any indirect, incidental, or consequential damages of any nature, or losses or expenses resulting from any defective product or the use of any product. Product(s) may not be resold, modified, or altered for resale without prior written approval from Eagle Biosciences, Inc.

For further information about this kit, its application or the procedures in this kit insert, please contact the Technical Service Team at Eagle Biosciences, Inc. at info@eaglebio.com or at 866-411-8023.



BIOMEDICA