

Bream GH Polyclonal Antibody (Rabbit)

Description

Bream growth hormone (GH) polyclonal antibody was prepared by immunising rabbits with recombinant bream GH (*Acanthopagrus butcheri*). The immunoglobulin fraction was purified by Protein A affinity chromatography. The antibody has 100% cross-reactivity with black bream GH but < 0.2% to trout GH, goldfish GH, human GH, rat GH and porcine GH. It is suitable for use in immunoassay, immunoblotting and immunohistochemistry

References:

Knibb W *et al* (1991) DNA Sequence **2**, 121-123

Immunogen:

Bream (*Acanthopagrus butcheri*) GH conjugated to ovalbumin

Source:

Rabbit antiserum. The immunoglobulin fraction was purified by Protein A affinity chromatography of serum obtained from a rabbit immunised with recombinant bream (*Acanthopagrus butcheri*) growth hormone.

Specificity:

Specific for Bream GH (*Acanthopagrus butcheri*). Cross-reactivity to trout GH (*Oncorhynchus mykiss*), goldfish GH (*Carassius auratus*), human GH, rat GH and porcine GH is < 0.2%.

Appearance:

Lyophilized powder

Storage/Stability:

At least 2 years at 2 - 4°C (lyophilized).
After reconstitution store at -20°C or -80°C
Avoid freeze-thaw cycles

Reconstitution:

Dissolve the vial contents in 200 µl sterile water or buffer solution and stand for 10 minutes at room temperature. Dilute further in a buffer suitable for the application.

Application and Titre:

Each application and titre should be determined in house.
As a guideline, following reconstitution a dilution in the range 1:500 – 1:32,000 should give a suitable absorbance in a conventional format.

Product Codes and Pricing

Bream GH
Polyclonal antibody

200 µg **PAG1**

Related Products:

Bream GH

****NOT FOR USE IN HUMANS****

GroPep Bioreagents Pty Ltd
51 West Thebarton Road
Thebarton SA 5031
Australia

ABN 93 147 032 166

Telephone: +61 8 7222 1051

Postal Address:
PO Box 10065
Adelaide Business Centre
SA 5000
Australia

Email:
info@gropep.com

Internet:
www.gropep.com

