



EAGLE
BIOSCIENCES

Anti-Human CYP450 3A4 IgG

Product Number: CYP30-A100

Aliquot: 100 µL

Storage: -80° C

PRODUCT SPECIFICATION SHEET

DESCRIPTION:	<u>This polyclonal antibody is produced in rabbit and is provided as antiserum.</u>
MOLECULAR WEIGHT:	~56 kDa
STORAGE BUFFER:	TBS
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS
IMMUNOGEN:	Peptide unique to human cytochrome P450 3A4
PURIFICATION:	None
CLASS:	Antiserum
APPLICATIONS:	Western blot: 1:2000 yields a positive result with 0.2 µg of recombinant cytochrome P450 3A4 using HRP conjugated anti-rabbit secondary antibody.
SPECIES SPECIFICITY:	Human.

Product Developed and Manufactured in the USA

EAGLE BIOSCIENCES, INC.
20A NW Blvd., Suite 112
Nashua, NH 03063

Phone: 617-419-2019
Fax: 617-419-1110
Email: info@EagleBio.com
www.EagleBio.com