



EAGLE
BIOSCIENCES

Anti-Human GST M1-1 Antisera

Catalog Number: GSM00-A100

Aliquot: 100 μ L

Storage: -80°C

PRODUCT SPECIFICATION SHEET

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum and is known not to cross-react with the pi or alpha class GST's.
MOLECULAR WEIGHT:	26.5 kDa
CONCENTRATION:	Not Applicable
STORAGE BUFFER:	Not Applicable
STORAGE:	-80°C. Aliquot once thawed and stored at -80°C. Avoid repeated freeze thaw cycles.
IMMUNOGEN:	Recombinant Human GST M1-1 expressed in <i>E. coli</i> .
PURIFICATION:	Not Purified
PURITY:	Not Purified
APPLICATIONS:	Western blot: 1:500 Is reactive with $\geq 0.2 \mu\text{g}$ of recombinant human GST M1-1 (product # GS 65) in conjunction with an HRP conjugated donkey anti-rabbit secondary antibody.
SPECIES SPECIFICITY:	Human and Rat
REFERENCE:	Holmes, W., et al. <i>Br. J. Cancer</i> , 62 : 209-212, (1990)

Product Developed and Manufactured in the USA

EAGLE BIOSCIENCES, INC.
20A NW Blvd., Suite 112
Nashua, NH 03063

Phone: 617-419-2019
Fax: 617-419-1110
Email: info@EagleBio.com
www.EagleBio.com