



**EAGLE**  
BIOSCIENCES

**Anti-Human Cytochrome  
P450 2D6**

Product Number: P4530-A100

Aliquot: 100  $\mu$ L

Storage: -20° C

**PRODUCT SPECIFICATION SHEET**

<b>Immunoglobulin Type:</b>	Polyclonal Rabbit antiserum
<b>Immunogen:</b>	Recombinant human P450 2D6 expressed in E. coli.
<b>Specificity:</b>	<u>WESTERN BLOTS:</u> Recognizes P450 2D6 in human liver microsomes. Reacts with $\geq 10 \mu\text{g}$ of recombinant human P450 2D6.
<b>Applications:</b>	Western blots.
<b>Quantity:</b>	100 $\mu$ L of rabbit antiserum. Sufficient for 10 Western blots.
<b>Working Dilution:</b>	A dilution of 1/1000 for Western blots will yield a positive result with $\geq 10 \mu\text{g}$ of recombinant human cytochrome P450 2D6 using alkaline phosphatase-conjugated anti-rabbit IgG secondary antibody according to the manufacturer's instructions.
<b>Storage:</b>	Store at -20°C. Do not freeze and thaw repeatedly. Aliquot the product for later use. Refreezing the product may result in loss of activity.

*Product Developed and Manufactured in the USA*

EAGLE BIOSCIENCES, INC.  
20A NW Blvd., Suite 112  
Nashua, NH 03063

Phone: 617-419-2019  
Fax: 617-419-1110  
Email: [info@EagleBio.com](mailto:info@EagleBio.com)  
[www.EagleBio.com](http://www.EagleBio.com)