



EAGLE
BIOSCIENCES

Anti-Human Prostaglandin E Synthase


Product Number: PES00-A100

Aliquot: 100 μ L

Storage: -80° C

DESCRIPTION:	Anti-Human Prostaglandin E Synthase
HOST:	Rabbit
STORAGE BUFFER:	DPBS
STORAGE:	-80°C, Avoid multiple freeze/thaw cycles
IMMUNOGEN:	Human Prostaglandin E Synthase
PURIFICATION:	Antiserum
APPLICATIONS:	For Western Blotting with ECL, use a 1:1000 dilution with the following conditions: Incubate with the primary antibody for 2 hours. Incubate with the secondary antibody for 1 hour. Extra washing is recommended to reduce non-specific background. Use an exposure time of 45 seconds.

PGES antibody test 30 sec exposure



* V = vector-transfected HCT-116 cells
P = PGES-transfected HCT-116 cells

Product Developed and Manufactured in the USA

EAGLE BIOSCIENCES, INC.
20A NW Blvd., Suite 112
Nashua, NH 03063

Phone: 617-419-2019
Fax: 617-419-1110
Email: info@EagleBio.com
www.EagleBio.com