



**EAGLE**  
BIOSCIENCES

## Anti-PGH Synthase

Catalog Number: PGS00-A100

Aliquot: 100  $\mu$ l

Storage: -70°C

### PRODUCT SPECIFICATION SHEET

<b>Immunoglobulin type:</b>	Polyclonal goat IgG.
<b>Specificity:</b>	Reacts with PGH Synthase from sheep and most other mammalian species including rodents.
<b>Uses:</b>	Immunoassay, Immunohistochemistry, and Western blots.
<b>Contents (per vial):</b>	100 $\mu$ L of goat antiserum containing 0.05% sodium azide.
<b>Stability:</b>	Stable for several weeks at 4° C, several months at -20° C, and indefinitely at -70° C. Avoid repeated freezing and thawing.
<b>Protein Concentration:</b>	Approximately 6.0 mg/mL.

#### References:

1. Smith, W.L. and Wilkin, G.P., Prostaglandins 13, 873-892 (1977).
2. Lee, Makau., Biochemica, 2, 19-21, (1996).

*Product Developed and Manufactured in the USA*

EAGLE BIOSCIENCES, INC.  
20A NW Blvd., Suite 112  
Nashua, NH 03063

Phone: 617-419-2019  
Fax: 617-419-1110  
Email: [info@EagleBio.com](mailto:info@EagleBio.com)  
[www.EagleBio.com](http://www.EagleBio.com)