



EAGLE
BIOSCIENCES

Anti-Protein Kinase C ϵ

Catalog Number: PKC00-A100

Aliquot: 100 μ L

Storage: -70°C

PRODUCT SPECIFICATION SHEET

| | |
|-----------------------------|--|
| DESCRIPTION: | This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C ϵ (PKC ϵ). |
| MOLECULAR WEIGHT: | 83.7 kDa Calculated |
| CONCENTRATION: | Not Determined |
| ISOTYPE: | NA |
| STORAGE BUFFER: | Phosphate Buffered Saline; 0.05% Sodium Azide |
| STORAGE: | -70°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month. |
| IMMUNOGEN: | KLH conjugated synthetic peptide corresponding to amino acids 721- 737 of PKC ϵ . Sequence: ⁷²¹ N-Q-E-E-F-K-G-F-S-Y-F-G-E-D-L-M-P ⁷³⁷ |
| PURIFICATION: | NA |
| APPLICATIONS: | Western Blotting: 1:1,000 Dilution |
| SPECIES SPECIFICITY: | Human, mouse and rat among other mammalian species. |
| REFERENCE: | JBC (1993) 268 :20110-20115 |

Product Developed and Manufactured in the USA

EAGLE BIOSCIENCES, INC.
20A NW Blvd., Suite 112
Nashua, NH 03063

Phone: 617-419-2019
Fax: 617-419-1110
Email: info@EagleBio.com
www.EagleBio.com